



Curie Point Injector JCI-22 VOC Analysis Cleanroom




mini-PAT

Quartz Glass Sample Tube
ID: 4mm Length: 36mm
Absorbent :Tenax GR 25mg
Pyrofoil :F315 (315c)



JCI-22
JAS-15M II

Curie Point Injector JCI-22 VOC in a car




JCI-22
JAS-15M II

Japan Analytical Industry Co., Ltd



Curie Point Injector JCI-22 VOC Analysis at Park

mini-PAT F280



JCI-22
JAS-15M II

Japan Analytical Industry Co., Ltd